

# Convenience with Versatility



Wilson Frame™



# Wilson Frame™



## Radiolucent Frame for Posterior Surgeries

The **Wilson Plus™** Radiolucent Frame provides a convenient and stable method of maintaining patients in a flexed position for Laminectomy, Decompression, Disc Surgery and Microdiscectomy procedures. The Wilson Plus Universal model is designed for use on a Trios or Modular Table System Imaging Top or on any general surgery table for ultimate versatility.

### Features

- Radiolucent design allows for unrestricted imaging at the surgical site
- Pads flex to adjust spinal lordosis and open the disc space for improved access
- Pads may be moved laterally to improve ventilation and relieve pressure on the abdomen
- ProneView® provides superior head and eye protection for the prone patient

### Specifications

<b>Weight of Frame</b>	
Model 5319G	40 lbs. (18 kg)
Model 5321G	50 lbs. (22 kg)
<b>Lateral Adjustment of Pads</b>	0–6 in. (0–15 cm)
<b>Intraoperative Flex of Pads for Lordosis Control</b>	2 in. (5 cm)
<b>Patient Weight Capacity</b>	300 lbs. (136 kg)

### Ordering Information

**REF**

**5319G Wilson Plus**  
Radiolucent Frame, Universal

**5321G Wilson Plus**  
Radiolucent Frame, Spinal Surgery Top

### Table Information

**REF**

7803 Trios Base  
7927 Trios Radiolucent Imaging Top

5803 Modular Table System Base  
5927 Radiolucent Imaging Top

### Individual Patient Care Kits

The patient care kits help protect the frame while promoting proper positioning and patient hygiene.

#### Kit Includes:

2 pad sleeves, 1 transverse bar pad, 1 base cover drape, 2 TLC™ Arm Cradles and either a ProneView® Cushion Insert Headrest Pillow or GentleTouch® Pillow.

**REF**

**5322 PV Wilson Plus**  
Patient Care Kit with ProneView (6/cs)


**5322 Wilson Plus**  
Patient Care Kit with GentleTouch (6/cs)



30031 Ahern Avenue  
Union City, CA 94587-1234 USA  
Telephone: 510-429-1500  
Toll Free: 800-777-4674  
Fax: 510-429-8500  
Outside USA: +1-510-429-1500  
[mizuho.com](http://mizuho.com)  
[newhipnews.com](http://newhipnews.com)

**EC REP**

Emergo Europe  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands



Note: Mizuho OSI is constantly improving its products. All specifications are subject to change without notice. • Mizuho OSI is a Delaware Corporation. Manufactured in the USA.

2016 ©Mizuho OSI **REF** NW0889 REV B