



new

ProneView®
Protective
Helmet System

Patient Positioning
Made Easier



ProneView[®] Protective Helmet System

Adjustable Protective Helmet – Makes Positioning Your Patients Easier

- Four large knobs are attached to the helmet that are easier to adjust once the patient has been rotated onto the surgical table
- Easier alignment on to the mirror platform
- Larger range of height adjustment for better neutrality of the neck
- Compatible with current disposable face cushions

The Mirror – Simple Yet Effective

- Easily aligns with the pegs located on the base of the Adjustable Protective Helmet
- Allows for continuous eye monitoring during procedure
- Fast and easy clean up, with nothing attached to the surface

Ordering Information

REF

D28502CE	ProneView Adjustable Protective Helmet
D28580CE	ProneView Mirror Platform with Post Holes
D28503CE	ProneView Large Cushion Inserts (10/cs)
D28507CE	ProneView Small Cushion Inserts (10/cs)

Note: For initial purchase, the helmet and mirror must both be ordered.
Replacement parts may be purchased separately.



Quick and Easy Adjustments

Offers a 1.65" (4.19 cm) range of elevation for easier patient positioning



Independent adjustability for improved neutrality of the patients neck



30031 Ahern Avenue
Union City, CA 94587-1234 USA
Telephone: 510-429-1500
Toll Free: 800-777-4674
Fax: 510-429-8500
Outside USA: +1-510-429-1500
www.mizuhoosi.com

Mizuho OSI's ISO13485 Quality Management System is certified by TÜV Rheinland[®] LGA Products GmbH.



Mizuho OSI is the exclusive worldwide distributor of the ProneView[®] Protective Helmet System. US Patents 6112333, 7426763 and 6490737. ProneView[®] is a registered trademark of Dupaco[®].