

STRETCHER AND WHEEL CHAIR TEMPUR-PEDIC® MEDICAL PADS

Tempur-Pedic® Medical Wheel Chair Pads
Pressure Equalization Pads, 3" multi-layer construction

REF	
69461	Tempur-Pedic® Medical Wheel Chair Pads 16"x16"x3" (41cm x 41 cm x 8 cm)
69462	Tempur-Pedic® Medical Wheel Chair Pads 18"x16"x3" (46 cm x 41 cm x 8 cm)
69463	Tempur-Pedic® Medical Wheel Chair Pads 18"x18"x3" (46 cm x 46 cm x 8 cm)
69464	Tempur-Pedic® Medical Wheel Chair Pads 20"x18"x3" (51cm x 46 cm x 8 cm)
69465	Tempur-Pedic® Medical Wheel Chair Pads 20"x16"x3" (51cm x 41 cm x 8 cm)

Tempur-Pedic® Medical Stretcher Pads

- Tempur-Pedic® Medical Stretcher Pads are 3" thick
- Optimal pressure equalization for patients up to 700 lbs (317 kg)
- Single soft Velcro® strip attaches to most stretcher models
- Angled and beveled corners improve access to transport handles and IV Pole receptacles

REF	
69024	Tempur-Pedic® Stretcher Pad 24" x 75.5"
69026	Tempur-Pedic® Stretcher Pad 26" x 75.5"
69265	Tempur-Pedic® Stretcher Pad 26.5" x 78"
69295	Tempur-Pedic® Stretcher Pad 29.5" x 75.5"

¹ DeFloor, T., DeBacquer, D., and Grypdonck, M. (2005). The Effect of Various Combinations of Turning and Pressure Reducing Devices on the Incidence of Pressure Ulcers, International Journal of Nursing Studies, 42 (2005), 37-46.

² Support Surfaces: Viscoelastic Expanded Urethane versus Conventional Foam
- Eric Flam, PhD, PE and Loretta Raab, RPh, CCP, ECPM Jan/Feb 2005, pp 29 - 34

³ Preventing Pressure Ulcers: An Evaluation of Four Operating-Table Mattresses Tom Defloor and Johan D. S. De Schuijmer
Copyright 2000 by W.B. Saunders Company

Mizuho OSI can provide Tempur-Pedic® Medical Pads to fit most surfaces. To provide your patients with the Pressure Equalization benefits of Tempur-Pedic® Medical Pads, please contact your local Mizuho OSI sales representative.




Distributed By
MIZUHO OSI
30031 AHERN AVENUE
UNION CITY, CA 94587/USA
BUSINESS: 510 - 429 -1500
INTERNATIONAL: 001- 510 - 429 - 1500
TOLL FREE: 800 - 777- 4674
FAX: 510 - 429 -8500
WWW.MIZUHOSI.COM
WWW.NEWHIPNEWS.COM

ISO certified by TÜV Rheinland® of North America



Mizuho OSI is the worldwide distributor of Tempur-Pedic® Medical Pressure Equalization Pads.
Note: Mizuho OSI is constantly improving its products. All specifications are subject to change without notice.
Mizuho OSI is a Delaware Corporation. Manufactured in the USA.
Tempur-Pedic® is a registered trademark of Tempur-Pedic® Medical, Inc.

 Tempur-Pedic North America, LLC • 1000 Tempur Way, Lexington, KY 40511 • Tel: 859-455-1000 • Int: 001-859-259-1000
Emergo Europe • Molenstraat 15 • 2513 BH The Hague • The Netherlands • Tel: (31) (0) 70 345-5870 • Fax: (31) (0) 70 346-7299



P R E S S U R E M A N A G E M E N T

TEMPUR-PEDIC® MEDICAL
PRESSURE EQUALIZATION PADS

TAKE THE PRESSURE OFF YOUR PATIENTS.®



TAKE THE PRESSURE OFF YOUR PATIENTS.®

IMPROVING PATIENT SAFETY AND COMFORT WITH PROGRESSIVE RESISTANCE TECHNOLOGY

HOW TEMPUR-PEDIC® MEDICAL PRESSURE EQUILIZATION PADS WORK

You can feel the difference in a Tempur-Pedic® Medical Pad just by touching it. It's called Progressive Resistance Technology. It means the more deeply a patient compresses a Tempur-Pedic® Medical Pad, the more it envelops the contact points with a fluid density that reduces peak pressures. This allows extraordinary weight support through pressure redistribution while remaining optimally comfortable to the patient.

IMPROVES PATIENT CARE

- Equalizes pressure load against the body
- Reduces incidence of nosocomial pressure ulcers¹
- Enhances patient comfort
- Adjusts firmness to body temperature
- Latex-free, anti-microbial materials
- Hypo-allergenic

SUPPORTS SURGICAL TEAM

- Material is radiolucent, MRI safe
- Waterfall edges
- Will not absorb or bleed energy
- Moisture resistant
- Non-skid bottom
- Easy cleaning with standard disinfectants
- Flame barrier complies with TB133

TEMPUR-PEDIC® MATERIAL VS CONVENTIONAL FOAM

Studies show Tempur-Pedic® material provides up to 4 times the pressure reduction of conventional polyurethane in controlled testing.²




Moisture impervious, vapor permeable membrane

Flame resistant barrier mattress assembly, CTB 133 certified

Tempur 110kg/cubic meter density for peak loading force distribution

Tempur-HD soft, dense absorption reduces localized pressure concentration



TEMPUR-PEDIC® MEDICAL PADS FOR MIZUHO OSI TABLES AND DEVICES

MTS Radiolucent Imaging Top

REF	
5840-7028	5927 Tempur-Pedic®Table Pad

MTS Spinal Surgery Top

REF	
6946	Pillow Standard Size (Sling Pad)
6947	Pillow Queen Size (Sling Pad)
5379-7722	5379 Articulating/5376 Translating, 1" Arm Board Pad
5840-7030	Support Section (Leg Board) Pad
5840-7569	SST Single Chest Pad

MTS Ortho Trauma Top

REF	
5855-7297	5855 Tempur-Pedic®Table Pad
5855-7367	Well Leg Tempur-Pedic®Pad
5855-7377	5855 Transfer Board Tempur-Pedic® Pad

hana®, hanaSSXT® and PROfx®

REF	
6875-7129	6875 Perineal Post Tempur-Pedic®Table Pad
6875-7120	6875 Main Table Tempur-Pedic®Table Pad
6850-7600	6850 PROfx®and hana® Arm Board Tempur-Pedic® Pad
6850-7590	6850 - 170 Leg Support Tempur-Pedic® Pad
6850-7560	6850 Perineal Post Tempur-Pedic®Table Pad
6850-7570	6850 Main Table Tempur-Pedic® Pad
6800-7140	6850-6000 PROfx®Auxiliary Imaging Top/6800 and 6862 Allegro®, Tempur-Pedic®Table Pad also fits Allegro® 6800 and 6862
6850-7580	6850-126 Lateral Perineal Support Tempur-Pedic®Pad
6870-39	hanaSSXT® Body Tempur-Pedic®Table Pad
6870-40	hanaSSXT® Perineal Tempur-Pedic®Table Pad

Allegro® Imaging Table

REF	
6800-7140	6800, 6862 Allegro® / 6850-6000 PROfx® Auxiliary Imaging Top, Tempur-Pedic® Table Pad
5363-7000	6900-12 Radiolucent Arm Board Tempur-Pedic®Pad
6900-25-7020	6900-25 AV Fistula Board Tempur-Pedic® Pad
6900-51-7001	6900-51 Universal Foot Board Tempur-Pedic® Pad

Other Replacement Pads for Discontinued Tables

REF	
6977-7000	Axis System Tempur-Pedic®Supine Imaging Top
6866-7700	6866 Tempur-Pedic®Table Pad
5895-7133	5895 Max Access™Tempur-Pedic®Table Pad
6800-7140	6862 Allegro Tempur-Pedic®Table Pad, also fits Allegro® 6800, PROfx® Auxiliary Imaging Top
Andrews Spinal TableTempur-Pedic® Pads available on request.	



TEMPUR-PEDIC® MEDICAL ACCESSORY PADS

Arm Board Pads

REF	
5356-7000	5356 Tempur-Pedic®Arm Board Pad, 2"
5366-7000	5356 Tempur-Pedic®Arm Board Pad, 3"
69444	5356 Tempur-Pedic®Arm Board Pad, 4"
5857-7000	5857 Cross Arm Support 1" Tempur-Pedic®Pad
5364-7018	Ultra Shoulder Translating Tempur-Pedic®Arm Board Pad

Miscellaneous Replacement and Positioning Pads

REF	
3369-1	Adult Perineal Post Tempur-Pedic®Perineal Post Pad 6" diameter
5806-7010	5806, Radiolucent Imaging Top Tempur-Pedic®Table Pad, 5'
5880-7150	5880-09 Supine Operating Top Tempur-Pedic®Pad
5847-218	5847/5848 Pelvic Reconstruction Kit Tempur-Pedic®Pad Set

Tempur-Pedic® Medical Pads and Positioners – 6pc set or individual

REF	
6950	6950-01 through 6950-06 Tempur-Pedic® O.R. Pad Set, 6 pc.

6950-01 Tempur-Pedic® Medical Knee Block A size 10" x 8" x 8"

REF

6950-02 Tempur-Pedic® Medical Knee Block B size 14" x 14"x 4"

REF

6950-03 Tempur-Pedic® Medical Basic Roll size 18"L x 6" diameter

REF

6950-04 Tempur-Pedic® Medical Spinal Wedge size 24" x 10" x 6.75"

REF

6950-05 Tempur-Pedic® Medical Cradle Block size 20" x 16" x 6"

REF

6950-06 Tempur-Pedic® Medical Rectangular Pad size 14" x 9" x 3"



TEMPUR-PEDIC® MEDICAL PADS FOR GENERAL SURGICAL TABLES AND DEVICES

Tempur-Pedic® Medical Pads have been shown to effectively equalize pressure forces over all contact surfaces with the patient. ³



All Tempur-Pedic® Medical General Surgical Table Pads are three piece sets with perineal cut-out. Individual pad sections can be ordered by placing an 'H' for Head, 'B' for Back or 'F' for Foot at the end of the pad set part number.

STERIS 3080, 3085 THREE PIECE PAD WITH PERINEAL CUT-OUT

REF	
69402	Three Piece Tempur-Pedic® Medical Pad w/Perineal Cut-Out, 2"
69403	Three Piece Tempur-Pedic® Medical Pad w/Perineal Cut-Out, 3"
694035	Three Piece Tempur-Pedic® Medical Pad w/Perineal Cut-Out, 3.5"

SKYTRON 6001, 6002, 6500 and 6600 THREE PIECE PAD WITH PERINEAL CUT-OUT

REF	
69412	Three Piece Tempur-Pedic® Medical Pad w/Perineal Cut-Out, 2"
69413	Three Piece Tempur-Pedic® Medical Pad w/Perineal Cut-Out, 3"
694135	Three Piece Tempur-Pedic® Medical Pad w/Perineal Cut-Out, 3.5"

SKYTRON 6700 and 6701 THREE PIECE PAD WITH PERINEAL CUT-OUT

REF	
69472	Three Piece Tempur-Pedic® Medical Pad w/Perineal Cut-Out, 2"
69473	Three Piece Tempur-Pedic® Medical Pad w/Perineal Cut-Out, 2"
694735	Three Piece Tempur-Pedic® Medical Pad w/Perineal Cut-Out, 3.5"
694935	

SKYTRON 3501 and 3600 THREE PIECE PAD WITH PERINEAL CUT-OUT

REF	
69482	Three Piece Tempur-Pedic® Medical Pad w/Perineal Cut-Out, 2"
69483	Three Piece Tempur-Pedic® Medical Pad w/Perineal Cut-Out, 3"
694935	Three Piece Tempur-Pedic® Medical Pad w/Perineal Cut-Out, 3.5"

SKYTRON 3500 and Amsco C-max THREE PIECE PAD WITH PERINEAL CUT-OUT

REF	
69492	Three Piece Tempur-Pedic® Medical Pad w/Perineal Cut-Out, 2"
69493	Three Piece Tempur-Pedic® Medical Pad w/Perineal Cut-Out, 3"
694935	Three Piece Tempur-Pedic® Medical Pad w/Perineal Cut-Out, 3.5"

BERCHTOLD B810 THREE PIECE PAD WITH PERINEAL CUT-OUT

REF	
69422	Three Piece Tempur-Pedic® Medical Pad w/Perineal Cut-Out, 2"
69423	Three Piece Tempur-Pedic® Medical Pad w/Perineal Cut-Out, 3"
694235	Three Piece Tempur-Pedic® Medical Pad w/Perineal Cut-Out, 3.5"

TEMPUR-PEDIC® MEDICAL REPLACEMENT PADS FOR ARM BOARDS

24" Flat Arm Board Pads

REF	
69441	1" x 24" Flat Arm Board Pad
69442	2" x 24" Flat Arm Board Pad
69443	3" x 24" Flat Arm Board Pad
69444	4" x 24" Flat Arm Board Pad

26" Flat Arm Board Pads

REF	
694461	1" x 26" Flat Arm Board Pad
694462	2" x 26" Flat Arm Board Pad
694463	3" x 26" Flat Arm Board Pad
694464	4" x 26" Flat Arm Board Pad

24" Contoured (Channel) Arm Board Pads

REF	
69452	2" x 24" Contoured Arm Board Pad
69453	3" x 24" Contoured Arm Board Pad
69454	4" x 24" Contoured Arm Board Pad

26" Contoured (Channel) Arm Board Pads

REF	
694562	2" x 26" Contoured Arm Board Pad
694563	3" x 26" Contoured Arm Board Pad
694564	4" x 26" Contoured Arm Board Pad

